



**Helio KIDS**  
 ORDER FORM - INTERNATIONAL

**1- Purchasing information**

P.O. number:\* \_\_\_\_\_ Date: \_\_\_\_\_  
 Company Name:\* \_\_\_\_\_  
 Purchasing contact:\* \_\_\_\_\_  
 Phone:\* \_\_\_\_\_  
 Email:\* \_\_\_\_\_  
 Contact preference:\* Phone  Email   
 Marked for: \_\_\_\_\_

**3- Shipping information**

Ship to:\* \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_  
 Zip Code: \_\_\_\_\_  
 Phone: \_\_\_\_\_

**2- Clinician information**

Clinician:\* \_\_\_\_\_  
 Phone:\* \_\_\_\_\_  
 Email:\* \_\_\_\_\_

Contact preference:\* Phone  Email   
 MA083  MA084  MA085  MA086  MA087  MA088  MA089   
*English manual (STD) German manual Italian manual Dutch manual Spanish manual French manual Portuguese manual*

MODEL			USD (taxes not included)
WC06	<input type="checkbox"/>	Ultralight Folding Carbon Fiber Helio KIDS Wheelchair Including fixed seat depth OR depth adjustable back (79,4 kg - 175 lb weight capacity)	\$4,395
WC10	<input type="checkbox"/>	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber	\$35

**FREE Growth Kit within the first 5 years of the Helio KIDS wheelchair ownership**  
 The growth kit includes carbon crossbraces, crossbrace links, seat sling, back upholstery, front riggings and footplates.

SEAT WIDTH <sup>1</sup>			NC
SW	<input type="checkbox"/>	30.48 (12") <input type="checkbox"/> 33.02 (13") <input type="checkbox"/> 35.56 (14") <input type="checkbox"/> 38.10 (15") <input type="checkbox"/> 40.64 (16") <input type="checkbox"/>	NC

[1. Width Calculator](#)

↓ **Caution: complete only one of the following two sections** ↓

SEAT DEPTH (Available with WC-19)			
SD12	<input type="checkbox"/>	30.48cm/12" Seat sling - frame of 35.56 cm/14"	NC
SD13	<input type="checkbox"/>	33.02.02cm/13" Seat sling - frame of 35.56 cm/14"	NC
SD14	<input type="checkbox"/>	35.56cm/14" Seat sling - frame of 35.56 cm/14"	NC
SD15	<input type="checkbox"/>	38.10cm/15" Seat sling - frame of 38.10cm/15"	NC
SD16	<input type="checkbox"/>	40.64cm/16" Seat sling - frame of 40.64cm/16"	NC

FIXED BACK CANE (MUST select ONE push handle options)			
BC11-1215	<input type="checkbox"/>	Straight back cane 30.48-38.10cm (12-15")	NC
BC11-1518	<input type="checkbox"/>	Straight back cane 38.10-45.72cm (15-18")	NC
BC11-1821	<input type="checkbox"/>	Straight back cane 45.72-53.34cm (18-21") STD	NC
BC12-1315	<input type="checkbox"/>	Back cane 8° bend 34.29-38.10cm (13.5-15")	NC
BC12-1518	<input type="checkbox"/>	Back cane 8° bend 38.10-45.72cm (15-18")	NC
BC12-1821	<input type="checkbox"/>	Back cane 8° bend 45.72-53.34cm (18-21") STD	NC

Specify height:<sup>5</sup> \_\_\_\_\_ measured from seat sling to top of back upholstery

ANGLE ADJUSTABLE BACK CANE (MUST select ONE push handle options)			
BC13-1518	<input type="checkbox"/>	Straight back cane 38.10-45.72cm (15-18")	NC
BC13-1821	<input type="checkbox"/>	Straight back cane 45.72-53.34cm (18-21") STD	NC
BC14-1518	<input type="checkbox"/>	Back cane 8° bend 38.10-45.72cm (15-18")	NC
BC14-1821	<input type="checkbox"/>	Back cane 8° bend 45.72-53.34cm (18-21") STD	NC

Select angle for the back cane: 85°  90°  95° STD  100°   
 105°  110°

Specify height:<sup>5</sup> \_\_\_\_\_


<sup>5-</sup> If not indicated, default height setting is at the lowest position

SEAT DEPTH- Depth and angle adustable (Not available with WC-19) <sup>3</sup>			
SD14	<input type="checkbox"/>	30.48 - 38.10 cm (12-15") (35.56/14" seat sling/frame)	NC
SD15	<input type="checkbox"/>	33.02 - 40.64 cm (13-16") (38.10/15" seat sling/frame)	NC
SD16	<input type="checkbox"/>	35.56 - 43.18 cm (14-17") (40.64/16" seat sling/frame)	NC
SD17	<input type="checkbox"/>	38.10 - 45.72 cm (15-18") (43.18/17" seat sling/frame)	NC
SD18	<input type="checkbox"/>	40.64 - 48.26 cm (16-19") (45.72/18" seat sling/frame)	NC

Specify depth<sup>4</sup>:

<sup>3</sup> - Maximum weight of 175 lb with a depth adjustable backrest

<sup>4</sup> - If not indicated, the depth will be equal to the seat sling



DEPTH AND ANGLE ADJUSTABLE BACK CANE (MUST select ONE push handle options)			
BC15-912	<input type="checkbox"/>	Straight back cane 22.86-30.48cm (9-12")	NC
BC15-1215	<input type="checkbox"/>	Straight back cane 30.48-38.10cm (12-15")	NC
BC15-1518	<input type="checkbox"/>	Straight back cane 38.10-45.72cm (15-18")	NC
BC15-1821	<input type="checkbox"/>	Straight back cane 45.72-53.34cm (18-21") STD	NC
BC16-1215	<input type="checkbox"/>	Back cane 8° bend 30.48-38.10cm (12-15")	NC
BC16-1518	<input type="checkbox"/>	Back cane 8° bend 38.10-45.72cm (15-18")	NC
BC16-1821	<input type="checkbox"/>	Back cane 8° bend 45.72-53.34cm (18-21") STD	NC

Select angle for the back cane: 90°  95° STD  100°  105°

Specify height:<sup>5</sup> \_\_\_\_\_

TRANSIT TIE DOWN			
AC04	<input type="checkbox"/>	Transit Tie Down (Unoccupied)	NC
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>6</sup> Short belt 50-104 cm (20" - 41")	\$365
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>6</sup> Long belt 50-116 cm (20" - 46")	\$365
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt <sup>6</sup>	\$335

<sup>6-</sup> Not available with depth adjustable back

## Helio Kids FRONT Seat-to-floor Height - Please select your choice

		7,62cm (3") caster			10,16cm (4") caster			12,7cm (5") caster			15,24cm (6") caster			17,78cm (7") caster			20,32cm (8") caster			
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			
Forks / Position		STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	
4" 10,2 cm	P1	33.02cm 13"	35.56cm 14"	38.10cm 15"	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"													
	P2	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"										
5" 12,7 cm	P1																			
	P2				35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"										
	P3	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"							
	P4	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"	40.64cm 16"	43.18cm 17"	45.72cm 18"	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"				
7" 17,8 cm	P1										41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"				
	P2							41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"				
	P3							43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	
	P4	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20"	53.3cm 21"	
	P5	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20"	53.3cm 21"	49.53cm 19.5"	52.07cm 20.5"	54.6cm 21.5"	

Height available +/- 0,65cm (1/4") (see stem bolt and fork combination at the end of the order form)

CASTERS			Availability				
CA01-03 to CA01-08	<input type="checkbox"/>	Composite wheel w/ PU tire (1")				3", 4", 5", 6", 7", 8"	NC
CA10 to CA11	<input type="checkbox"/>	Pneumatic casters (1 1/4")				6", 8"	\$175
CA14 to CA16	<input type="checkbox"/>	Newton UltraCaster Composite (1 1/2")				4", 5", 6"	\$98
CA17 to CA19	<input type="checkbox"/>	Newton UltraCaster Aluminum (1 1/2")				4", 5", 6"	\$198
CASTER OPTIONS							
CA12	<input type="checkbox"/>	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) <sup>1</sup>					\$320
CA13	<input type="checkbox"/>	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)					\$515

1 - Not available with WC-19

### KIDS REAR Seat-to-floor Height

WHEEL SIZES		Height available		Specify desired height: _____
<input type="checkbox"/>	20" (451) wheel	31.75cm / 12.5"	to 43.81cm / 17.25"	
<input type="checkbox"/>	22" (501) wheel	34.29cm / 13.5"	to 46.35cm / 18.25"	
<input type="checkbox"/>	24" (540) wheel	36.83cm / 14.5"	to 48.89cm / 19.25"	

Note:

Footrest						
SWING-IN/SWING-OUT FOOTREST						
FR03-01	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 34.29 - 43.18cm (13.5" - 17")				Set at: * NC
FR03-02	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5")				Set at: * NC
FR03-03	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 39.29 - 52.07cm (15.5" - 20.5")				Set at: * NC
FR04-01	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 34.29 - 43.18cm (13.5" - 17")				Set at: * NC
FR04-02	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5")				Set at: * NC
FR04-03	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 39.37 - 52.07cm (15.5" - 20.5")				Set at: * NC
FR15-01	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 34.29 - 40.64cm (13.5" - 16") <sup>5</sup>				Set at: * NC
FR15-02	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5") <sup>5</sup>				Set at: * NC
FR15-03	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 39.37 - 52.07cm (15.5" - 20.5") <sup>5</sup>				Set at: * NC
FR06	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Standard (no extension) <sup>1</sup>		19.05 - 27.94cm (7.5" - 11")		Set at: * NC
FR07	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Short extension 2" <sup>1,2</sup>		24.13 - 33.02cm (9.5" - 13")		Set at: * NC
FR08	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Medium extension 4" <sup>1,2</sup>		24.13 - 36.83cm (9.5" - 14.5")		Set at: * NC
FOOTREST OPTIONS						
FR09	<input type="checkbox"/>	High mount (footplate attached to the hanger with clamps) <sup>3</sup>		12.70 - 34.29cm (5" - 13.5")		Set at: * \$85
FR10	<input type="checkbox"/>	Elevating legrest with calf pad, length adjustable <sup>4</sup> Left <input type="checkbox"/> Right <input type="checkbox"/>		35.56 - 55.88 cm (14"-22")		Set at: * \$248
FR99	<input type="checkbox"/>	Omit Footrest				NC

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount and the elevating legrest

\* +/- 1/2" (1.27 cm)

4- The elevating legrest replace the footrest side selected; Not compatible with seat width of 12" and 13" and one piece footplate; 5- Not available with 8" caster

FOOTPLATES							
FP01	<input type="checkbox"/>	Standard - Composite <sup>1</sup>					NC
FP02	<input type="checkbox"/>	Newton adjustable angle - Composite					\$230
FP03	<input type="checkbox"/>	Newton adjustable angle - Aluminum <sup>2</sup>					\$255
FP04	<input type="checkbox"/>	Newton adjustable angle - Carbon fiber with rubber inserts					\$450
FP23	<input type="checkbox"/>	Newton adjustable angle extra deep - Aluminum					\$235
FP06	<input type="checkbox"/>	One piece adjustable angle flip-up - Aluminum					\$205
FP07	<input type="checkbox"/>	One piece adjustable angle flip-up - Carbon fiber with rubber inserts					\$470
FOOTPLATE ACCESSORIES							
FP08	<input type="checkbox"/>	Length adjustable heel loop (not available with one piece footplate)					\$45
FP09	<input type="checkbox"/>	Calf strap					\$40
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap					\$125

1- Not available with front mount; 2- The Newton adjustable angle (Aluminum) is the only available option with a 12 and 13" wide frame

SWING AWAY RESIDUAL LIMB SUPPORT							
FR27	<input type="checkbox"/>	Residual limb support - flat 15.24x20.32cm (6"x8")	Left <input type="checkbox"/>	Right <input type="checkbox"/>			\$440
FR28	<input type="checkbox"/>	Residual limb support - flat 20.32x24.4cm (8"x10")	Left <input type="checkbox"/>	Right <input type="checkbox"/>			\$440
FR29	<input type="checkbox"/>	Residual limb support - flat 25.4x25.4cm (10"x10")	Left <input type="checkbox"/>	Right <input type="checkbox"/>			\$440
FR30	<input type="checkbox"/>	Residual limb support - flat 25.4x35.56cm (10"x14")	Left <input type="checkbox"/>	Right <input type="checkbox"/>			\$440
FR32	<input type="checkbox"/>	Residual limb support - contour 15.24x20.32cm (6"x8")	Left <input type="checkbox"/>	Right <input type="checkbox"/>			\$440
FR33	<input type="checkbox"/>	Residual limb support - contour 20.32x24.4cm (8"x10")	Left <input type="checkbox"/>	Right <input type="checkbox"/>			\$440

PUSH HANDLES OPTIONS (ONE option at a time)							
PH01	<input type="checkbox"/>	Push handles					STD
PH05	<input type="checkbox"/>	Fold down push handles <sup>1</sup>					\$140
PH06-920	<input type="checkbox"/>	clamp on height adjustable push handles 22.86-50.80cm (9-20") <sup>1,2</sup>					\$330
PH06-923	<input type="checkbox"/>	clamp on height adjustable push handles 22.86-58.40cm (9-23") <sup>1,2</sup>					\$330
PH02	<input type="checkbox"/>	Omit push handles					NC

STROLLER HANDLE AND STABILIZER BAR (only available with push handle)							
PH04	<input type="checkbox"/>	Stroller handle					\$175
PH03	<input type="checkbox"/>	Newton folding stabilizer bar (not available with stroller handle)					\$275

1- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock; 2- Not compatible with depth adjustable back, 8" bend back canes and amputee axle plates

BACK UPHOLSTERY							
BH07	<input type="checkbox"/>	Slip-on Active-Mesh back upholstery					NC
BH08	<input type="checkbox"/>	Adjustable tension Active-Mesh back upholstery					\$180
BH99	<input type="checkbox"/>	Omit Back upholstery					STD

SEAT CUSHION							
SC01	<input type="checkbox"/>	5.08 cm (2") Cushion					\$180
SC02	<input type="checkbox"/>	7.62 cm (3") Cushion					\$200
SC99	<input type="checkbox"/>	Omit seat cushion					STD

### Rear wheels

ONE ARM DRIVE <sup>1</sup> Only available with Newton One rear wheels, permanent axles, 0° camber							
HR22	<input type="checkbox"/>	One arm drive - aluminum inner handrim <sup>2</sup>	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		\$1,510
HR60	<input type="checkbox"/>	One arm drive - plastic coated inner handrim <sup>2</sup>	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		\$1,595
HR98	<input type="checkbox"/>	Omit handrim on opposite side of the One Arm Drive (OAD)					NC

REAR WHEELS				Available									
RW01	<input type="checkbox"/>	Mag <sup>3</sup>				20", 22" 24"	\$57						
RW02	<input type="checkbox"/>	Newton One - spoke wheel <sup>4</sup>				20", 22" 24"	NC						
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel				20", 22", 24"	\$335						
RW04	<input type="checkbox"/>	Spinergy Spox	Blk STD <input type="checkbox"/>	Wht <sup>5</sup> <input type="checkbox"/>	Blue <sup>5</sup> <input type="checkbox"/>	Red <sup>5</sup> <input type="checkbox"/>	Yel <sup>5</sup> <input type="checkbox"/>	22", 24"	\$865				
RW05	<input type="checkbox"/>	Spinergy LX	Blk STD <input type="checkbox"/>	Pur <sup>5</sup> <input type="checkbox"/>	Wht <sup>5</sup> <input type="checkbox"/>	Blue <sup>5</sup> <input type="checkbox"/>	Red <sup>5</sup> <input type="checkbox"/>	Yel <sup>5</sup> <input type="checkbox"/>	Pink <sup>5</sup> <input type="checkbox"/>	Grn <sup>5</sup> <input type="checkbox"/>	Org <sup>5</sup> <input type="checkbox"/>	22", 24"	\$1,095
RW07	<input type="checkbox"/>	Spinergy CLX <sup>6</sup>	Blk STD <input type="checkbox"/>	Pur <sup>5</sup> <input type="checkbox"/>	Wht <sup>5</sup> <input type="checkbox"/>	Blue <sup>5</sup> <input type="checkbox"/>	Red <sup>5</sup> <input type="checkbox"/>	Yel <sup>5</sup> <input type="checkbox"/>	Pink <sup>5</sup> <input type="checkbox"/>	Grn <sup>5</sup> <input type="checkbox"/>	Org <sup>5</sup> <input type="checkbox"/>	24"	\$1,600

TIRES				Available			
<b>Solid</b>							
TI01	<input type="checkbox"/>	Soft Urethane 3.49cm (1 3/8")				20", 22" 24"	STD
TI02	<input type="checkbox"/>	Soft Urethane knobby 3.49cm (1 3/8 ") <sup>7</sup>				22", 24"	\$105
TI04	<input type="checkbox"/>	Soft Urethane 2.54cm (1") - Shox				22", 24"	\$125
<b>Pneumatic</b>							
TI05	<input type="checkbox"/>	Pneumatic 3.49cm (1 3/8")				20", 22", 24"	NC
TI08	<input type="checkbox"/>	Pneumatic with Airless insert 3.49cm (1 3/8") <sup>8</sup>				22", 24"	\$90
TI06	<input type="checkbox"/>	Pneumatic HP, Low Tread, Puncture resistant, black 2.54cm (1") - SpeedLite				20", 22", 24"	\$110
TI07	<input type="checkbox"/>	Pneumatic HP, Medium Tread, Puncture resistant, black 2.54cm (1") - TrailBlazer				24"	\$140

1- Offered with handrims on each side of the wheelchair; Select outer handrim in the handrims section; 3- Not compatible with HP tires; Only available with aluminum or plastic coated handrim, standard position;

4- Not compatible with Shox tires; 5- Colour spokes require an additional 10 days lead-time; 6- Only available with 1" tires; 7- The 22" is only available with Mags and Newton One rear wheels; 8- Not available with Spinergy CLX

## Handrims

HANDRIM			Available		
HR01	<input type="checkbox"/>	Aluminum anodized	20", 22", 24"		STD
HR61	<input type="checkbox"/>	Aluminum black hard anodized <sup>2</sup>	20", 22", 24"		\$75
HR02	<input type="checkbox"/>	Plastic coated <sup>1</sup>	20", 22", 24"		\$110
HR03	<input type="checkbox"/>	High friction coated <sup>1,2,3</sup>	22", 24"		\$380
HR07-HR18 and HR65-HR71	<input type="checkbox"/>	Newton Air Grip <sup>2</sup> <span style="color: red;">See on page 6 for color choice and price</span>	20", 22", 24"		
HR28	<input type="checkbox"/>	Surge <sup>1,2,3</sup> (1 1/2" Oval handrim w/ Gripton high friction strip)	22", 24"		\$425
HR05	<input type="checkbox"/>	Surge LT <sup>1,2,3</sup> (1 1/8" Oval handrim w/ Gripton high friction strip)	22", 24"		\$465
HR04	<input type="checkbox"/>	Natural Fit <sup>1,2,3</sup>	22", 24"		\$335
HR62	<input type="checkbox"/>	Simi H <sup>1,2,3</sup> (1 1/4" Oval handrim with high friction strip)	22", 24"		\$395
HR63	<input type="checkbox"/>	Nova H <sup>1,2,3</sup> (1 1/4" Oval handrim with high friction strip)	22", 24"		\$425
HR64	<input type="checkbox"/>	Optimum H <sup>1,2,3</sup> (1 1/8" heart-shaped handrim)	22", 24"		\$325
HR99	<input type="checkbox"/>	Omit handrim			NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;

## Handrims option

ASSEMBLY POSITION					
HR29	<input type="checkbox"/>	Standard Handrim position <span style="float: right;">Full tab</span>			STD
HR030	<input type="checkbox"/>	Narrow Handrim position no cut <span style="float: right;">Full tab</span>			NC
HR30	<input type="checkbox"/>	Narrow Handrim position with cut <sup>3</sup> <span style="float: right;">Short tab 2 positions</span>			\$32
HR31	<input type="checkbox"/>	Super Narrow Handrim <sup>3</sup> <span style="float: right;">Short tab 2 positions STD <input type="checkbox"/> Short tab 1 position <input type="checkbox"/></span>			\$32

THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS					
HR23	<input type="checkbox"/>	Standard grip			\$95
HR27	<input type="checkbox"/>	Super grip			\$150
HR099	<input type="checkbox"/>	Omit thumb piece			NC


3- Only available with Aluminum and Air Grip handrims;

AXLE					
AX01	<input type="checkbox"/>	Quick release axle <sup>1</sup>			STD
AX02	<input type="checkbox"/>	Permanent axle			NC
AX03	<input type="checkbox"/>	Quad release axle <sup>1</sup>			\$135

AXLE PLATE					
AX04	<input type="checkbox"/>	Amputee adjustable axle plate <sup>2</sup>			\$210

1- Not required with one arm drive 2- May interfere with swing away armrest, depth adjustable back and limit seat to floor height available

REAR WHEEL CAMBER					
CM00	<input type="checkbox"/>	0° camber			STD
CM03	<input type="checkbox"/>	3° camber			NC
CM06	<input type="checkbox"/>	6° camber (not available with amputee axle plate)			\$165

CENTER OF GRAVITY																																																																						
	<p><b>Choose position of rear wheel axle plate:</b></p> <p>Measure distance X from front of back post to center of axle receiver.</p>																																																																					
	X =	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 10%;">Most stable</th> <th style="width: 10%;">3,2</th> <th style="width: 10%;">3,8</th> <th style="width: 10%;">4,4</th> <th style="width: 10%;">5,1</th> <th style="width: 10%;">5,7</th> <th style="width: 10%;">6,3</th> <th style="width: 10%;">7</th> <th style="width: 10%;">7,6</th> <th style="width: 10%;">8,3</th> <th style="width: 10%;">8,9</th> <th style="width: 10%;">9,5</th> <th style="width: 10%;">10,1</th> <th style="width: 10%;">Most responsive</th> </tr> </thead> <tbody> <tr> <td>AMPUTEE<sup>1</sup></td> <td>2,5</td> <td>3,2</td> <td>3,8</td> <td>4,4</td> <td>5,1</td> <td>5,7</td> <td>6,3</td> <td>7</td> <td>7,6</td> <td>8,3</td> <td>8,9</td> <td>9,5</td> <td>10,1</td> </tr> <tr> <td></td> <td>1"</td> <td>1/4"</td> <td>1/2"</td> <td>3/4"</td> <td>2"</td> <td>1/4"</td> <td>1/2"</td> <td>3/4"</td> <td>3"</td> <td>1/4"</td> <td>1/2"</td> <td>3/4"</td> <td>4"</td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </tbody> </table>	Most stable	3,2	3,8	4,4	5,1	5,7	6,3	7	7,6	8,3	8,9	9,5	10,1	Most responsive	AMPUTEE <sup>1</sup>	2,5	3,2	3,8	4,4	5,1	5,7	6,3	7	7,6	8,3	8,9	9,5	10,1		1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>												
	Most stable	3,2	3,8	4,4	5,1	5,7	6,3	7	7,6	8,3	8,9	9,5	10,1	Most responsive																																																								
	AMPUTEE <sup>1</sup>	2,5	3,2	3,8	4,4	5,1	5,7	6,3	7	7,6	8,3	8,9	9,5	10,1																																																								
	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"																																																									
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																									
<p><b>If not indicated, default setting is at 3,8 cm (1.5")</b></p> <p>Possible conflict between a seat depth under 40,6 cm (16") and a center of gravity over 6,3 cm (2.5").</p> <p>Possible conflict between a 66 cm (26") rear wheel and a center of gravity over 6,3 cm (2.5").</p> <p>Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG</p> <p>1-Only available with the amputee axle plate</p>																																																																						

WHEEL LOCK					
WL08	<input type="checkbox"/>	Newton Wheel lock - push to lock			STD
WL14	<input type="checkbox"/>	6" removable extension handle for push to lock			\$83
WL16	<input type="checkbox"/>	Aluminum push to lock with extensions			\$83
WL02	<input type="checkbox"/>	Pull to lock			NC
WL04	<input type="checkbox"/>	Newton grade aid push to lock (works best with pneumatic or high pressure tire)			\$165
WL13	<input type="checkbox"/>	6" removable extension handle for pull to lock and grade aid			\$83
WL05	<input type="checkbox"/>	Scissors Lock			\$155
WL06	<input type="checkbox"/>	Attendant Lock (additional set of wheel lock mounted to the back)			\$215
WL15	<input type="checkbox"/>	Unilateral wheel lock - pull to lock <sup>1</sup> <span style="float: right;">Left <input type="checkbox"/> Right <input type="checkbox"/></span>			\$335
WL99	<input type="checkbox"/>	Omit wheel locks			NC

1- For more information about incompatibilities and restriction, refer to the table at page 6 and the "How-to documents" page of our website:


[https://www.motioncomposites.com/en\\_intl/support-and-education/how-to-documents](https://www.motioncomposites.com/en_intl/support-and-education/how-to-documents)

Armrest										
<b>HEIGHT ADJUSTABLE "T" ARMREST<sup>1</sup></b>										
AR13-610	<input type="checkbox"/>	25.4cm (10") Desk length armpad <sup>4</sup>	15.25 - 25.40cm (6" - 10")						\$350	
AR03-812	<input type="checkbox"/>	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")						\$350	
AR03-1014	<input type="checkbox"/>	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")						\$350	
AR14-610	<input type="checkbox"/>	35.56cm (14") Full length armpad <sup>4</sup>	15.25 - 25.40cm (6" - 10")						\$350	
AR04-812	<input type="checkbox"/>	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")						\$350	
AR04-1014	<input type="checkbox"/>	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")						\$350	
<b>LOCKING FLIP-UP ARMREST<sup>2</sup></b>										
AR07	<input type="checkbox"/>	Locking Flip-Up Armrest								\$355
<b>ARMREST OPTIONS</b>										
AR12	<input type="checkbox"/>	35.56cm (14") Full Length Armpad - Gel Ovation (pair) <sup>3</sup>								\$165
AR99	<input type="checkbox"/>	Omit armrest								NC
<b>1- May have 1/2" adjustment restriction on both limits of the range; 2- May interfere with rigid back hardware; Only available with 18"-21" straight aluminum back canes; Not available with back upholstery</b>										
<b>3- Not available with locking flip-up armrest and short "T" armrest; 4- May have some restriction with seat-to-floor height and rear wheel</b>										
<b>SIDE GUARDS (not available with "T" Armrest)</b>										
SG01	<input type="checkbox"/>	Plastic side guards (removable)								\$210
SG03	<input type="checkbox"/>	Carbon fiber side guards (removable)								\$315
SG06	<input type="checkbox"/>	Plastic fender side guards - straight bracket (removable)								\$295
SG08	<input type="checkbox"/>	Plastic fender side guards - offset bracket (removable) <sup>1</sup>								\$295
SG04	<input type="checkbox"/>	Carbon fiber fender side guards (removable)								\$550
SG99	<input type="checkbox"/>	Omit side guard								NC
<b>SIDE GUARDS OPTION</b>										
SG07	<input type="checkbox"/>	Side guards clip (pair)								NC
<b>1- only compatible with 0° camber; Increase distance between fender side guards by 1 1/8"</b>										
<b>POSITIONING BELT</b>										
SB01-S	<input type="checkbox"/>	Velcro adjustable - Standard 88.9cm (35")								\$50
SB01-L	<input type="checkbox"/>	Velcro adjustable - Long 152.40cm (60")								\$50
SB02-S	<input type="checkbox"/>	Auto buckle - Standard 88.9cm (35")								\$95
SB02-L	<input type="checkbox"/>	Auto buckle - Long 152.40cm (60")								\$95
SB04-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Medium 142.24cm (56")								\$95
SB04-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Long 154.94cm (61")								\$95
SB05-S	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Small 129.54cm (51")								\$185
SB05-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Medium 142.24cm (56")								\$185
SB05-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Long 154.94cm (61")								\$185
SB06-S	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Small 129.54cm (51")								\$205
SB06-M	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Medium 142.24cm (56")								\$205
SB06-L	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Long 154.94cm (61")								\$205
SB03	<input type="checkbox"/>	Hardware attachment for 2 extra points (1 pair)								\$17
SB08-M	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 73.66 (18")								\$310
SB08-L	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 93.98cm (24")								\$310
SB99	<input type="checkbox"/>	Omit positioning belt								STD
<b>HEADREST<sup>1</sup></b>										
AC01-S	<input type="checkbox"/>	Profiled Headrest - Small 8.8x15.2cm (3.5"x6")								\$195
AC01-M	<input type="checkbox"/>	Profiled Headrest - Medium 11.4x17.8cm (4.5"x7")								\$195
AC01-L	<input type="checkbox"/>	Profiled Headrest - Large 14x2.3cm(5.5"x8")								\$195
AC28	<input type="checkbox"/>	Headrest mounting hardware with removable non folding stabilizer bar <sup>1</sup>								\$320
<b>1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar</b>										
<b>ACCESSORIES</b>										
AC03	<input type="checkbox"/>	Anti-tippers								\$150
AC29	<input type="checkbox"/>	Swing-away anti-tipper <sup>1</sup>	Left <input type="checkbox"/>	Right <input type="checkbox"/>					\$165	
AC27	<input type="checkbox"/>	Tip assist <sup>2</sup>								\$67
AC06	<input type="checkbox"/>	Cane and crutch holder	Left <input type="checkbox"/>	Right <input type="checkbox"/>					\$130	
AC07	<input type="checkbox"/>	Oxygen tank holder <sup>4</sup> - D and E Type - installed on the left								\$335
AC08	<input type="checkbox"/>	Spoke guards <sup>3</sup>	Available: 20", 22", 24", 26"							\$125
AC09	<input type="checkbox"/>	Footrest neoprene impact guard (pair) - Black								\$100
AC10	<input type="checkbox"/>	Seat pouch								\$83
AC30	<input type="checkbox"/>	Newton multi-purpose pouch								\$62
AC23	<input type="checkbox"/>	Newton Backpack								\$130
AC26	<input type="checkbox"/>	BodyPoint Mobility Bag								\$135
AC12	<input type="checkbox"/>	Velcro Closing Strap								\$57
AC13	<input type="checkbox"/>	Tool kit								\$50
AC16	<input type="checkbox"/>	Motion Composites Polo Shirt	W/S <input type="checkbox"/>	W/M <input type="checkbox"/>	W/L <input type="checkbox"/>				\$36	
			M/S <input type="checkbox"/>	M/M <input type="checkbox"/>	M/L <input type="checkbox"/>	M/XL <input type="checkbox"/>	M/XXL <input type="checkbox"/>	M/XXXL <input type="checkbox"/>		
<b>1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 34,3cm (13 1/2") and 50,2cm (19 3/4"); 2- Cannot be combined with Anti Tippers</b>										
<b>3- Not available with Mag wheels and Spineray; 4- Not compatible with swing away armrest or amputee axle plate and WC-19 transit tie-down</b>										

Colors																					
Options <sup>1</sup>	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
AirGrip	<input type="checkbox"/> \$110		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	\$265
Frame <sup>2</sup>		<input type="checkbox"/> STD	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	NC
FC99	<input type="checkbox"/>	Custom Color	<a href="#">Specify PPG Paintit color:</a>														<i>additional 10 days lead-time</i>	\$390			

- 1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;  
 2- with black hydroformed aluminum crossbrace

### Fork and Stem bolt combinations



### Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
22 in	17	1.0 to 3.25
	18+	1.0 to 3.75
	16	1.0 to 2.0
24 in	17	1.0 to 2.75
	18+	1.0 to 3.75

Notes :