



VELOCE
 ORDER FORM - INTERNATIONAL

1- Purchasing information

P.O. number:* _____ Date: _____
 Company Name:* _____
 Purchasing contact:* _____
 Phone:* _____
 Email:* _____
 Contact preference:* Phone Email

Marked for: _____

3- Shipping information

Ship to:* _____
 Name: _____
 Address: _____
 City: _____ State: _____
 Zip Code: _____
 Phone: _____

2- ATP/Clinician information

ATP/Clinician:* _____
 Phone:* _____
 Email:* _____
 Contact preference:* Phone Email * Required


MA083 <input type="checkbox"/>	MA084 <input type="checkbox"/>	MA085 <input type="checkbox"/>	MA086 <input type="checkbox"/>	MA087 <input type="checkbox"/>	MA088 <input type="checkbox"/>	MA089 <input type="checkbox"/>
English Manual (STD)	German Manual	Italian Manual	Dutch Manual	Spanish Manual	French Manual	Portuguese Manual

MODEL								USD (Taxes not included)						
WC20	<input type="checkbox"/>	Ultralight VELOCE Wheelchair (120 kg - 265 lb weight limit)							\$5,195					
WC09	<input type="checkbox"/>	HD Heavy Duty package (159 kg - 350 lb weight capacity)							\$540					
WC10	<input type="checkbox"/>	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber							\$35					
Frame														
SEAT WIDTH ³ - cm (inches)														
SW	<input type="checkbox"/>	35.56 (14")	<input type="checkbox"/>	38.10 (15")	<input type="checkbox"/>	40.64 (16")	<input type="checkbox"/>	43.18 (17")	<input type="checkbox"/>	45.72 (18") ²	<input type="checkbox"/>	48.26 (19") ²	<input type="checkbox"/>	NC
	<input type="checkbox"/>	50.80 (20") ²				<input type="checkbox"/>	53.34 (21") ⁷	<input type="checkbox"/>	55.88 (22") ⁷		<input type="checkbox"/>			\$535
SEAT DEPTH - cm (inches)														
SD	<input type="checkbox"/>	35.56 (14") ¹	<input type="checkbox"/>	38.10 (15")	<input type="checkbox"/>	40.64 (16")	<input type="checkbox"/>	43.18 (17")	<input type="checkbox"/>	45.72 (18")	<input type="checkbox"/>	48.26 (19")	<input type="checkbox"/>	NC
	<input type="checkbox"/>	50.80 (20")				<input type="checkbox"/>								\$515
FRAME LENGTH (Seat sling depth plus distance from seat sling to front of frame) ⁴														
SP01	<input type="checkbox"/>	Standard frame length configuration ⁵											STD	
SP02	<input type="checkbox"/>	Frame length 2.54 cm (1") longer than standard configuration											NC	
SP03	<input type="checkbox"/>	Frame length 5.08 cm (2") longer than standard configuration											NC	
SP04	<input type="checkbox"/>	Frame length 2.54 cm (1") shorter than standard configuration											NC	
SP05	<input type="checkbox"/>	Frame length 5.08 cm (2") shorter than standard configuration											NC	
TRANSIT TIE DOWN														
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt ⁶					Short belt 50-104 cm (20" - 41")				\$365			
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt ⁶					Long belt 50-116 cm (20" - 46")				\$365			
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt ⁶											\$335	

1- Only available with frame length 1" and 2" longer; 2- Available with HD option [3. Width calculator](#) 4- See chart on page 7
 5- 4.375" from seat sling to front frame; 6- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option; 7- Only available with HD option the one piece adjustable angle flip-up

FRONT Seat-to-floor height

Caster sizes	Height available +/- 1/4"	Specify desired height:
<input type="checkbox"/> 3" Caster	35.56cm to 53.34cm (14" to 21")	_____ (inches)
<input type="checkbox"/> 4" Caster	36.83cm to 53.34cm (14.5" to 21")	
<input type="checkbox"/> 5" Caster	39.37cm to 53.34cm (15.5" to 21")	
<input type="checkbox"/> 6" Caster	41.91cm to 53.34cm (16.5" to 21") (N/A with 90° front frame)	



CASTERS					Availability	
CA01-03 to CA01-06	<input type="checkbox"/>	Composite wheel w/ PU tire (1")			3", 4", 5", 6"	NC
CA10	<input type="checkbox"/>	Pneumatic casters (1 ¼")			6"	\$175
CA14 to CA16	<input type="checkbox"/>	Newton UltraCaster Composite (1 ½ ")			4", 5", 6"	\$98
CA17 to CA19	<input type="checkbox"/>	Newton UltraCaster Aluminum (1 ½ ")			4", 5", 6"	\$198
CASTER OPTIONS						
CA13	<input type="checkbox"/>	Carbon fiber Frog Legs suspension forks <i>(Not available with caster pin locks)</i>				\$515


REAR Seat-to-floor Height

WHEEL SIZES		Height available	Specify desired height: _____ (inches)
<input type="checkbox"/>	22" (501) wheel	34.29cm to 46.35cm (13.5" to 18.25")	
<input type="checkbox"/>	24" (540) wheel	36.83cm to 48.89cm (14.5" to 19.25")	
<input type="checkbox"/>	25" (559) wheel	38.10cm to 50.16cm (15" to 19.75")	
<input type="checkbox"/>	26" (590) wheel	39.37cm to 51.43cm (15.5" to 20.25")	


Configure

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1,25 cm (½") of selected Seat to Floor height"


FRONT FRAME ANGLE					
FA70	<input type="checkbox"/>	70° Front frame angle			NC
FA75	<input type="checkbox"/>	75° Front frame angle			NC
FA80	<input type="checkbox"/>	80° Front frame angle			NC
FA85	<input type="checkbox"/>	85° Front frame angle			NC
FA90	<input type="checkbox"/>	90° Front frame angle			NC



SEAT TO FOOTREST LENGTH			
Set Length to ¹ :		Length must be at least 4,5 cm (1 ¾") shorter than front seat-to-floor height	
1 - +/- 1.27 cm (1/2")			



Footrest width / Taper selection							
Seat Width		Straight ¹		1" Taper			
Outside Measurement		6,4 cm Narrower than Seat Width		11,4 cm Narrower than Seat Width			
35.56	14"	29.21 cm	11.5"	24.13 cm	9.5"		
38.10	15"	31.75 cm	12.5"	26.67 cm	10.5"		
40.64	16"	34.29 cm	13.5"	29.21 cm	11.5"		
43.18	17"	36.83 cm	14.5"	31.75 cm	12.5"		
45.72	18"	39.37 cm	15.5"	34.29 cm	13.5"		
48.26	19"	41.91 cm	16.5"	36.83 cm	14.5"		
50.80	20"	44.45 cm	17.5"	39.37 cm	15.5"		
53.34	21"	46.99 cm	18.5"	41.91 cm	16.5"		
55.88	22"	49.53 cm	19.5"	44.45 cm	17.5"		



¹ - Not compatible with auto-folding footplate

FOOTPLATES						
FP02	<input type="checkbox"/>	Newton adjustable angle footplate - Composites <i>(Two piece)</i>				\$230
FP03	<input type="checkbox"/>	Newton adjustable angle - Aluminum <i>(Two piece)</i>				\$255
FP04	<input type="checkbox"/>	Newton adjustable angle footplate - Carbon fiber with rubber inserts <i>(Two piece)</i>				\$450
FP06	<input type="checkbox"/>	One piece adjustable angle flip-up - Aluminum				NC
FP07	<input type="checkbox"/>	One piece adjustable angle flip-up - Carbon fiber with rubber inserts				\$570
FP28	<input type="checkbox"/>	One piece auto-folding adjustable angle - Aluminum				\$310
FP29	<input type="checkbox"/>	One piece auto-folding adjustable angle - Carbon fiber with rubber inserts				\$380
FOOTPLATE ACCESSORIES						
FP08	<input type="checkbox"/>	Length adjustable heel loop <i>(not available with one piece footplate)</i>				\$45
FP09	<input type="checkbox"/>	Calf strap				\$40
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap				\$125

Backrest

FIXED HEIGHT WITHOUT PUSH HANDLE OPTIONS						
BC03-9		Ultralight carbon fiber back cane	22.86cm (9")	*Temporarily unavailable		NC
BC03-12		Ultralight carbon fiber back cane	30.48cm (12")	*Temporarily unavailable		NC
BC03-15		Ultralight carbon fiber back cane	38.10cm (15")	*Temporarily unavailable		NC
BC03-18		Ultralight carbon fiber back cane	45.72cm (18")	*Temporarily unavailable		NC
BC03-21		Ultralight carbon fiber back cane	53.34cm (21")	*Temporarily unavailable		NC
BC04-13		8° Bend Carbon fiber back cane	34.29cm (13 ½)	*Temporarily unavailable		\$60
BC04-15		8° Bend Carbon fiber back cane	38.10cm (15")	*Temporarily unavailable		\$60
BC04-18		8° Bend Carbon fiber back cane	45.72cm (18")	*Temporarily unavailable		\$60

*** Replacement options below: 1-Carbon fiber height adjustable; 2- Fixed height aluminum; 3- carbon fiber height and angle adjustable**

HEIGHT ADJUSTABLE WITH PUSH HANDLES (MUST select ONE push handle options)						
BC06-912	<input type="checkbox"/>	Straight Carbon fiber back cane	22.86 - 30.48cm (9-12")	increment of 1/2"	Set at ² :	NC
BC06-1215	<input type="checkbox"/>	Straight Carbon fiber back cane	30.48 - 38.10cm (12-15")	increment of 1/2"	Set at ² :	NC
BC06-1518	<input type="checkbox"/>	Straight Carbon fiber back cane	38.10 - 45.72cm (15-18")	increment of 1/2"	Set at ² :	NC
BC06-1821	<input type="checkbox"/>	Straight Carbon fiber back cane	45.72 - 53.37cm (18-21")	increment of 1/2"	Set at ² :	NC
BC08-1315	<input type="checkbox"/>	8° Bend Carbon fiber back cane	34.29 - 38.10cm (13 ½-15")	increment of 1/2"	Set at ² :	\$60
BC08-1518	<input type="checkbox"/>	8° Bend Carbon fiber back cane	38.10 - 45.72 (15-18")	increment of 1/2"	Set at ² :	\$60
BC08-1821	<input type="checkbox"/>	8° Bend Carbon fiber back cane	45.72 - 53.34cm (18-21")	increment of 1/2"	Set at ² :	\$60

FIXED HEIGHT WITHOUT PUSH HANDLE OPTIONS						
BC30	<input type="checkbox"/>	Straight aluminum back cane	9"			NC
BC36	<input type="checkbox"/>	Straight aluminum back cane	12"			NC
BC42	<input type="checkbox"/>	Straight aluminum back cane	15"			NC
BC63	<input type="checkbox"/>	8° bend aluminum back cane	13 ½"			\$60
BC64	<input type="checkbox"/>	8° bend aluminum back cane	15"			\$60

ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)						
BC09-1215	<input type="checkbox"/>	Straight Carbon fiber back cane	30.48 - 38.10cm (12-15")	increment of 1/2"	Set at ² :	\$220
BC09-1518	<input type="checkbox"/>	Straight Carbon fiber back cane	38.10 - 45.72cm (15-18")	increment of 1/2"	Set at ² :	\$220
BC09-1821	<input type="checkbox"/>	Straight Carbon fiber back cane	45.72 - 53.37cm (18-21")	increment of 1/2"	Set at ² :	\$220
BC10-1315	<input type="checkbox"/>	8° bend Carbon fiber back cane	34.29 - 38.10cm (13 ½-15")	increment of 1/2"	Set at ² :	\$275
BC10-1518	<input type="checkbox"/>	8° bend Carbon fiber back cane	38.10 - 45.72 (15-18")	increment of 1/2"	Set at ² :	\$275
BC10-1821	<input type="checkbox"/>	8° bend Carbon fiber back cane	45.72 - 53.34cm (18-21")	increment of 1/2"	Set at ² :	\$275
		Select back cane angle:	85° ¹ <input type="checkbox"/>	90° <input type="checkbox"/>	95° STD <input type="checkbox"/>	100° <input type="checkbox"/>
					105° <input type="checkbox"/>	110° <input type="checkbox"/>

PUSH HANDLES OPTIONS (ONE option at a time)						
PH01	<input type="checkbox"/>	Push handles				NC
PH05	<input type="checkbox"/>	Fold down push handles ^{3,5}				\$140
PH06-920	<input type="checkbox"/>	Clamp on height adjustable push handles 22.86 - 50.80cm (9" - 20") ^{3,4,5}				\$330
PH06-923	<input type="checkbox"/>	Clamp on height adjustable push handles 22.86 - 58.4cm (9" - 23") ^{3,4,5}				\$330
PH02	<input type="checkbox"/>	Omit push handles				NC

1- Not available with a 8° bend back cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 120 kg/265lb weight limit

4- Not compatible with swing away armrests, 8° bend back canes; 5- Not compatible with headrest and attendant lock

BACK UPHOLSTERY						
BH07	<input type="checkbox"/>	Slip-on Active-Mesh back upholstery				NC
BH08	<input type="checkbox"/>	Adjustable tension Active-Mesh back upholstery				\$180
BH99	<input type="checkbox"/>	Omit Back upholstery				NC

SEAT CUSHION						
SC01	<input type="checkbox"/>	5.08 cm (2") Cushion				\$180
SC02	<input type="checkbox"/>	7.62 cm (3") Cushion				\$200
SC99	<input type="checkbox"/>	Omit seat cushion				STD

Rear wheels

ONE ARM DRIVE ¹		Only available with Newton One rear wheels, permanent axles, 0° camber											
HR22	<input type="checkbox"/>	One arm drive - aluminum inner handrim ^{2,3}	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>				\$1,510				
HR60	<input type="checkbox"/>	One arm drive - plastic coated inner handrim ^{2,3}	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>				\$1,595				
HR98	<input type="checkbox"/>	Omit handrim on opposite side of the One Arm Drive (OAD)							NC				
REAR WHEELS									Available				
RW01	<input type="checkbox"/>	Mag ⁴					22" 24"		\$57				
RW02	<input type="checkbox"/>	Newton One - spoke wheel ⁵					22" 24"		NC				
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel					22", 24", 25", 26"		\$335				
RW04	<input type="checkbox"/>	Spinergy Spox	Blk STD <input type="checkbox"/>	Wht ⁶ <input type="checkbox"/>	Blue ⁶ <input type="checkbox"/>	Red ⁶ <input type="checkbox"/>	Yel ⁶ <input type="checkbox"/>	22", 24", 25", 26"	\$865				
RW05	<input type="checkbox"/>	Spinergy LX	Blk STD <input type="checkbox"/>	Pur ⁶ <input type="checkbox"/>	Wht ⁶ <input type="checkbox"/>	Blue ⁶ <input type="checkbox"/>	Red ⁶ <input type="checkbox"/>	Yel ⁶ <input type="checkbox"/>	Pink ⁶ <input type="checkbox"/>	Grn ⁶ <input type="checkbox"/>	Org ⁶ <input type="checkbox"/>	22", 24", 25", 26"	\$1,095
RW07	<input type="checkbox"/>	Spinergy CLX ⁷	Blk STD <input type="checkbox"/>	Pur ⁶ <input type="checkbox"/>	Wht ⁶ <input type="checkbox"/>	Blue ⁶ <input type="checkbox"/>	Red ⁶ <input type="checkbox"/>	Yel ⁶ <input type="checkbox"/>	Pink ⁶ <input type="checkbox"/>	Grn ⁶ <input type="checkbox"/>	Org ⁶ <input type="checkbox"/>	24", 25"	\$1,600
TIRES									Available				
Solid													
TI01	<input type="checkbox"/>	Soft Urethane 3.49cm (1 3/8")					22" 24"		NC				
TI02	<input type="checkbox"/>	Soft Urethane knobby 3.49cm (1 3/8 ") ⁸					22", 24", 25", 26"		\$105				
TI04	<input type="checkbox"/>	Soft Urethane 2.54cm (1") - <i>Shox</i>					22", 24", 25", 26"		\$125				
Pneumatic													
TI05	<input type="checkbox"/>	Pneumatic 3.49cm (1 3/8")					22", 24", 25", 26"		NC				
TI08	<input type="checkbox"/>	Pneumatic with Airless insert 3.49cm (1 3/8") ⁹					22", 24", 26"		\$90				
TI06	<input type="checkbox"/>	Pneumatic HP, Low Tread, Puncture resistant, black 2.54cm (1") - <i>SpeedLite</i>					22", 24", 25", 26"		\$110				
TI07	<input type="checkbox"/>	Pneumatic HP, Medium Tread, Puncture resistant, black 2.54cm (1") - <i>TrailBlazer</i>					24", 25", 26"		\$140				

- 1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; 2- Only standard handrim position available; 3- Not available with HD option
 4- Not compatible with HP tires; Only available with aluminum or plastic coated handrim, standard position; 5- Not compatible with Shox tires; 6- color spokes require an additional 10 days lead-time;
 7- Only available with 1" tires; 8- 22" is only available with Mags and Newton One rear wheels; 9- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels

Handrims

HANDRIM								Available		
HR01	<input type="checkbox"/>	Aluminum anodized					22", 24", 25", 26"		STD	
HR61	<input type="checkbox"/>	Aluminum black hard anodized ²					22", 24", 25"		\$75	
HR02	<input type="checkbox"/>	Plastic coated ¹					22", 24", 25", 26"		\$110	
HR03	<input type="checkbox"/>	High friction coated ^{1,2,3}					22", 24", 25", 26"		\$380	
HR07-HR18 and HR65-HR71	<input type="checkbox"/>	Newton Air Grip ²	<i>See page 6 for the color selection and price</i>					22", 24", 25", 26"		
HR28	<input type="checkbox"/>	Surge ^{1,2,3} (1 1/2" Oval handrim w/ Gripton high friction strip)					22", 24", 25", 26"		\$425	
HR05	<input type="checkbox"/>	Surge LT ^{1,2,3} (1 1/8" Oval handrim w/ Gripton high friction strip)					22", 24", 25", 26"		\$465	
HR04	<input type="checkbox"/>	Natural Fit ^{1,2,3}					22", 24", 25", 26"		\$335	
HR62	<input type="checkbox"/>	Simi H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)					22", 24", 25", 26"		\$395	
HR63	<input type="checkbox"/>	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)					22", 24", 25", 26"		\$425	
HR64	<input type="checkbox"/>	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)					22", 24", 25", 26"		\$325	
HR99	<input type="checkbox"/>	Omit handrim							NC	

- 1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;

Handrims option

ASSEMBLY POSITION									
HR29	<input type="checkbox"/>	Standard Handrim position	Full tab						STD
HR030	<input type="checkbox"/>	Narrow Handrim position no cut	Full tab						NC
HR30	<input type="checkbox"/>	Narrow Handrim position with cut ³	Short tab 2 positions						\$32
HR31	<input type="checkbox"/>	Super Narrow Handrim ³	Short tab 2 positions STD <input type="checkbox"/>	Short tab 1 position <input type="checkbox"/>					\$32
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS									
HR23	<input type="checkbox"/>	Standard grip							\$95
HR27	<input type="checkbox"/>	Super grip							\$150
HR099	<input type="checkbox"/>	Omit thumb piece							NC

- 3- Only available with Aluminum and Air Grip handrims;


AXLE ¹

AX01	<input type="checkbox"/>	Quick release axle							STD
AX03	<input type="checkbox"/>	Quad release axle							\$135

- 1- Not required with one arm drive

REAR WHEEL CAMBER				
CM00	<input type="checkbox"/>	0° camber		STD
CM03	<input type="checkbox"/>	3° camber		NC
CM06	<input type="checkbox"/>	6° camber		\$165

CENTER OF GRAVITY



Choose position of rear wheel axle plate:
Measure distance X from front of back post to center of axle receiver.

Most stable												
2,5	3,2	3,8	4,4	5,1	5,7	6,3	7	7,6	8,3	8,9	9,5	10,2
1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

If not indicated, default setting is at 3,8 cm/1.5"
Possible conflict between a seat depth under 40.6 cm (16") and a center of gravity over 6,3 cm (2.5").
Possible conflict between a 26" rear wheel and a center of gravity over 6,3 cm (2.5").
Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

WHEEL LOCK				
WL08	<input type="checkbox"/>	Newton Wheel lock - push to lock		STD
WL14	<input type="checkbox"/>	6" removable extension handle for push to lock		\$83
WL16	<input type="checkbox"/>	Aluminum push to lock with extensions		\$83
WL02	<input type="checkbox"/>	Pull to lock		NC
WL04	<input type="checkbox"/>	Newton grade aid push to lock (works best with pneumatic or high pressure tire)		\$165
WL13	<input type="checkbox"/>	6" removable extension handle for pull to lock and grade aid		\$83
WL05	<input type="checkbox"/>	Scissors Lock		\$155
WL15	<input type="checkbox"/>	Unilateral wheel lock - pull to lock ¹	Left <input type="checkbox"/> Right <input type="checkbox"/>	\$335
WL99	<input type="checkbox"/>	Omit wheel locks		NC

1- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_intl/support-and-education/how-to-documents

Handrims				
HEIGHT ADJUSTABLE FLIP BACK ARMREST (convertible to "T" armrest) ¹				
AR01-812	<input type="checkbox"/>	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")	Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/> \$380
AR01-1014	<input type="checkbox"/>	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")	Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/> \$380
AR02-812	<input type="checkbox"/>	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")	Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/> \$380
AR02-1014	<input type="checkbox"/>	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")	Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/> \$380
HEIGHT ADJUSTABLE "T" ARMREST ¹				
AR13-610	<input type="checkbox"/>	25.4cm (10") Desk length armpad ⁴	15.25 - 25cm (6" - 10")	\$350
AR03-812	<input type="checkbox"/>	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")	\$350
AR03-1014	<input type="checkbox"/>	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")	\$350
AR14-610	<input type="checkbox"/>	35.56cm (14") Full length armpad ⁴	15.25 - 25cm (6" - 10")	\$350
AR04-812	<input type="checkbox"/>	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")	\$350
AR04-1014	<input type="checkbox"/>	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")	\$350
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²				
AR05	<input type="checkbox"/>	Tubular armrest - Short pad	20.32 - 27.94cm (8 - 11")	\$240
AR23	<input type="checkbox"/>	Tubular armrest - Long pad	20.32 - 27.94cm (8 - 11")	\$240
AR06	<input type="checkbox"/>	25.4cm (10") Desk Length flat armpad	20.32 - 27.94cm (8 - 11")	\$240
AR17	<input type="checkbox"/>	35.56cm (14") Full Length flat armpad	20.32 - 27.94cm (8 - 11")	\$240
ARMREST OPTIONS				
AR12	<input type="checkbox"/>	35.56cm (14") Full Length Armpad - Gel Ovation (pair) ³		\$165
AR99	<input type="checkbox"/>	Omit armrest		NC

1- May have a 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid back rests; Not available with fixed carbon fiber side guards

3- Only available with full length armrest; Not available with tubular (short and long) and short armrest; 4- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS (not available with Flip Back and "T" Armrest)				
SG01	<input type="checkbox"/>	Plastic side guards (removable)		\$210
SG02	<input type="checkbox"/>	Carbon fiber side guards (fixed)		\$275
SG03	<input type="checkbox"/>	Carbon fiber side guards (removable)		\$315
SG06	<input type="checkbox"/>	Plastic fender side guards - straight bracket (removable)		\$295
SG08	<input type="checkbox"/>	Plastic fender side guards - offset bracket (removable) ¹		\$295
SG04	<input type="checkbox"/>	Carbon fiber fender side guards (removable)		\$550
SG99	<input type="checkbox"/>	Omit side guard		NC
SIDE GUARDS OPTION				
SG07	<input type="checkbox"/>	Side guards clip (pair) ²		NC

1- only compatible with 0° camber; Increase distance between fender side guards by 1 1/8"; 2- Not available with fixed carbon fiber side guard

POSITIONING BELT			
SB01-S	<input type="checkbox"/>	Velcro adjustable - Standard 88.9cm (35") ¹	\$50
SB01-L	<input type="checkbox"/>	Velcro adjustable - Long 152.40cm (60")	\$50
SB02-S	<input type="checkbox"/>	Auto buckle - Standard 88.9cm (35") ¹	\$95
SB02-L	<input type="checkbox"/>	Auto buckle - Long 152.40cm (60")	\$95
SB04-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Medium 142.24cm (56")	\$95
SB04-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Long 154.94cm (61")	\$95
SB05-S	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Small 129.54cm (51")	\$185
SB05-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Medium 142.24cm (56")	\$185
SB05-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Long 154.94cm (61")	\$185
SB06-S	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Small 129.54cm (51")	\$205
SB06-M	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Medium 142.24cm (56")	\$205
SB06-L	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Long 154.94cm (61")	\$205
SB03	<input type="checkbox"/>	Hardware attachment for 2 extra points (1 pair)	\$17
SB08-M	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 73.66 (18")	\$310
SB08-L	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 93.98cm (24")	\$310
SB99	<input type="checkbox"/>	Omit positioning belt	STD

1- Not available with wheelchair width of 45.72cm (18") and more

HEADREST ¹			
AC01-S	<input type="checkbox"/>	Profiled Headrest - Small 8.8x15.2cm (3.5"x6")	\$195
AC01-M	<input type="checkbox"/>	Profiled Headrest - Medium 11.4x17.8cm (4.5"x7")	\$195
AC01-L	<input type="checkbox"/>	Profiled Headrest - Large 14x2.3cm(5.5"x8")	\$195
AC28	<input type="checkbox"/>	Headrest mounting hardware with tension bar	\$320

1- Only available with push handles

ACCESSORIES			
AC03	<input type="checkbox"/>	Anti-tippers	\$150
AC29	<input type="checkbox"/>	Swing-away anti-tipper ¹ Left <input type="checkbox"/> Right <input type="checkbox"/>	\$165
AC27	<input type="checkbox"/>	Tip assist ²	\$67
AC06	<input type="checkbox"/>	Cane and crutch holder Left Right <input type="checkbox"/>	\$130
AC08	<input type="checkbox"/>	Spoke guards ³ Available: 22", 24", 26"	\$125
AC10	<input type="checkbox"/>	Seat pouch	\$83
AC30	<input type="checkbox"/>	Newton multi-purpose pouch	\$62
AC26	<input type="checkbox"/>	BodyPoint Mobility Bag	\$135
AC23	<input type="checkbox"/>	Newton Backpack	\$130
AC12	<input type="checkbox"/>	Velcro Closing Strap	\$57
AC13	<input type="checkbox"/>	Tool kit	\$50
AC16	<input type="checkbox"/>	Motion Composites Polo Shirt W/S <input type="checkbox"/> W/M <input type="checkbox"/> W/L <input type="checkbox"/> M/S <input type="checkbox"/> M/M <input type="checkbox"/> M/L <input type="checkbox"/> M/XL <input type="checkbox"/> M/XXL <input type="checkbox"/> M/XXXL <input type="checkbox"/>	\$36

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 34,3cm (13 1/2") and 50,2cm (19 3/4"); 2- Cannot be combined with Anti Tippers

3- Not available with Mag wheels and Spinergy

Colors																					
Options ¹																					
AirGrip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	\$265
Frame	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	NC
Highlight ²	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	NC
FC99	<input type="checkbox"/>	Custom Color	Specify PPG Paintit color:																	additional 10 days lead-time	\$390

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- To see different highlight, visit our website: https://www.motioncomposites.com/en_intl/products/folding-wheelchairs/veloce

Frame length chart

Frame length	Seat Depth						
	35.56(14")	38.10(15")	40.64(16")	43.18(17")	45.72(18")	48.26(19")	50.80(20")
STD	N/A	38.10(15")	40.64(16")	43.18(17")	45.72(18")	48.26(19")	50.80(20")
1	38.10(15")	40.64(16")	43.18(17")	45.72(18")	48.26(19")	50.80(20")	N/A
2	40.64(16")	43.18(17")	45.72(18")	48.26(19")	50.80(20")	N/A	N/A
-1	N/A	N/A	38.10(15")	40.64(16")	43.18(17")	45.72(18")	48.26(19")
-2	N/A	N/A	N/A	38.10(15")	40.64(16")	43.18(17")	45.72(18")

Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.75 to 4
22 in	16+	2.5 to 4
24 in	16+	1 to 4

Notes :